



T: 0161 998 9219

Fast, Effective, Affordable Service

With many years' experience in the manufacturing and technical support of mass spectrometers, our engineering staff bring a wide range of expertise offering a fast and efficient service and repair of instruments. Official UK Distributor for Scion and Varian Instruments.

Exciting show for Mass Spec at the 2016 BMSS Meeting



The annual BMSS meeting held last month was once again a great success. This year the location was sunny Eastbourne, the event showcased new sessions and three excellent speakers.

The meeting also featured a full Exhibition, involving over 30 industrial groups and Mass Spec once again exhibited at the exhibition. This year we took the SCION 436-GC and SCION 456-GC instruments, as Mass Spec Technology Solutions is the UK distributor of SCION and GCMS instruments. SCION GCs provide the ultimate capability for today's analytical laboratory. The ultimate in GC platforms can be configured to handle virtually any GC analysis from the most basic single injector/detector to the most complex multi-valve/multi-column. A new platform designed to deliver top performance today and inspire further developments tomorrow. This incorporates legacy instruments of Bruker, Varian, GC and GC/MS.

[Read More](#)



Innovative Lab Bench NGenMS20

Mass Spec also exhibited their innovative Lab Bench, **NGenMS20**. The innovative design provides a vibration free bench while containing both vacuum pumps and compressed nitrogen generator thus saving lab space.

The new **NGenMS20** bench provides a compact solution to lab space and noise suppression requirements. An integrated CMS based nitrogen generator and vacuum pump noise suppression enclosures save valuable lab floor space and solves lab background noise issues. An 8 socket, anti-surge power distribution panel and twin lockable castors is supplied as standard and the bench has a high specification noise insulated backing pump enclosure.

[Read More](#)



Peter Heaps, Managing Director at Mass Spec comments “the conference was fantastic and we are already looking forward to next year’s event which will be held in Manchester I believe so perfect for us!! Thanks to everyone who dropped by our stand for a chat”

Mass Spec Technology Solutions provides total service support for mass spectrometry laboratories. Our experienced service team is committed to fast, effective service at affordable cost with a variety of plans to suit individual laboratory requirements.

[Click here to view our services](#)

Visit www.mass-spec.com for more information

T: 0161 998 9219 E: sales@mass-spec.com

Unit 16 Bentinck St, Ashton-Under-Lyne, Lancashire, OL7 0PT


Official UK Distributor

SERVICE • MAINTENANCE • TECH SUPPORT • SPARES • CONSUMABLES • TRAINING • MANUFACTURE



Copyright © 2016 Mass Spec Technology Solutions Ltd, All rights reserved.

Want to change how you receive these emails?
You can [update your preferences](#) or [unsubscribe from this list](#)

